



## LIFTING HANDLE

- Allows you to carry 4" and 6" specimens with one hand.
- Remove 6" and 4" specimens from curing tank.
- Safe and reduces back stress.

**CAT# TC6LH:** For 6" specimen.

**CAT# TC4LH:** For 4" specimen.



## 6CH ECONO CARRY HANDLE

The handy 6CH Econo Cylinder Carrier will fit any 6" (152mm) diameter plastic cylinder mold with a wide lip and can be used with lids in place.

Flexible plastic ring slides onto mold from bottom prior to use and is held in place by top lip. When ready to transport, carrier handles are rotated up to snap together in the middle.

- Carriers are packaged in cartons of 50.
- Recommended for single use.

**CAT# ECH**

**Made in the USA.**

